


Summer 2021 Trainer Catalogue



sales@buchananorthotics.co.uk 0141 440 1999 • 603 Helen Street, Glasgow G51 3AR

Coast













Dale



Shown in Black Mesh and Light Beige Suede

Arrow



Shown in Camel softee with Camel Suede and Tan detailing

Glen



Shown in Silver Mesh and Grey Nubuck





Take a step in the right direction...













Featur	res: round	toe,	2 x	3mm	poron	insoles	&	footb

	Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
		\leftrightarrow		\$	2		
5	1	57	85	217-235	201-219	24-33	223
Standard	2	58	88	224-242	208-226	25-34	231
=	3	60	90	230-248	214-232	26-35	239
žž.	4	61	93	237-255	221-239	26-35	248
0,	5	63	95	243-261	227-245	27-36	257
	6	64	98	250-268	234-252	28-37	265
	7	66	100	256-274	240-258	29-38	274
	8	67	103	263-281	247-265	29-38	282
	9	69	105	269-287	253-271	30-39	291
	10	70	108	276-294	260-278	31-40	300

Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
	\leftrightarrow	Θ				
1	61	91	229-247	209-227	24-33	223
2	62	94	236-254	216-234	25-34	231
3	64	96	242-260	222-240	26-35	239
4	65	99	249-267	229-247	26-35	248
5	67	101	255-273	235-253	27-36	257
6	68	104	262-280	242-260	28-37	265
7	70	106	268-286	248-266	29-38	274
8	71	109	275-293	255-273	29-38	282
9	73	111	281-299	261-279	30-39	291
10	74	114	288-306	268-286	31-40	300

Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
	\leftrightarrow	\bigcirc	4	4		
1	53	78	205-223	193-212	24-33	223
2	54	82	212-230	200-218	25-34	231
3	56	84	218-236	206-224	26-35	239
4	57	87	225-243	213-231	26-35	248
5	59	89	231-249	219-237	27-36	257
6	60	92	238-256	226-244	28-37	265
7	62	94	244-262	232-250	29-38	274
8	63	97	251-269	239-257	29-38	282
9	65	99	257-275	245-263	30-39	291
10	66	102	264-282	252-270	31-40	300

	Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
		\leftrightarrow		\$	2		
	1	57	85	229-247	201-209	30-39	223
0	2	58	88	236-254	208-226	31-40	231
eep	3	60	90	242-260	214-232	32-41	239
	4	61	93	249-267	221-239	32-41	248
	5	63	95	255-270	227-245	33-42	257
	6	64	98	262-280	234-252	34-43	265
	7	66	100	268-286	240-258	35-44	274
	8	67	103	275-293	247-265	35-44	282
	9	69	105	281-299	253-271	36-45	291
	10	70	108	288-306	260-278	37-46	300

Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
	Θ	Θ	A	2		
1	61	91	241-259	209-227	30-39	223
2	62	94	248-266	216-234	31-40	231
3	64	96	254-272	222-240	32-41	239
4	65	99	261-279	229-247	32-41	248
5	67	101	267-285	235-253	33-42	257
6	68	104	274-292	242-260	34-43	265
7	70	106	280-298	248-266	35-44	274
8	71	109	287-305	255-273	35-44	282
9	73	111	293-311	261-279	36-45	291
10	74	114	300-318	268-286	37-46	300
10	74	114	300-318	268-286	37-46	

Donna (Narrower Heel)

- (1	(Natiowel ricel)							
Si		Heel Width	Instep Circum.					
		Θ	2					
	1	-	-					
	2	54	198-217					
	2 3 4	56	204-223					
		57	211-230					
	5	59	217-236					
	6	60	224-243					
	7	62	230-249					
- 1	В	63	237-256					

Diana	
(Shallower	Forepart)

		Joint Circum.	Toe Depth
1		\$	
1	1	-	-
	2	214-233	21-30
]	3	220-239	22-31
4	4	227-246	22-31
4	5	233-252	23-32
4	6	240-259	24-33
4	7	246-265	25-34
J	8	253-272	25-34

	Heel Width	Instep Circum.	Size	Joint Circum.	Toe Depth
	Θ	2		\$	
1	-	-	1	-	-
2	58	206-225	2	226-245	21-30
3	60	212-231	3	232-251	22-31
4	61	219-238	4	239-258	22-31
5	63	225-244	5	245-264	23-32
6	64	232-251	6	252-271	24-33
7	66	238-257	7	258-277	25-34
8	67	245-264	8	265-284	25-34

		Heel Width	Instep Circum.	Size	Joint Circum.	Toe Depth
		\leftrightarrow	2		4	
ĺ	1	-	-	1	-	-
ĺ	2	50	190-209	2	202-221	21-30
ĺ	3	52	196-215	3	208-227	22-31
l	4	53	203-222	4	215-234	22-31
Į	5	55	209-228	5	221-240	23-32
l	6	56	216-235	6	228-247	24-33
ļ	7	58	222-241	7	234-253	25-34
	8	59	229-248	8	241-260	25-34

	Heel Width	Instep Circum.	Size	Joint Circum.	Toe Depth
	Θ	2		4	
1	-	-	1	-	-
2	50	190-209	2	202-221	21-30
3	52	196-215	3	208-227	22-31
4	53	203-222	4	215-234	22-31
5	55	209-228	5	221-240	23-32
6	56	216-235	6	228-247	24-33
7	58	222-241	7	234-253	25-34
8	59	229-248	8	241-260	25-34

tep Circum.	Size	Joint Circum.	Toe Depth
2		4	
-	1	-	-
190-209	2	202-221	21-30
196-215	3	208-227	22-31
203-222	4	215-234	22-31
209-228	5	221-240	23-32
216-235	6	228-247	24-33
222-241	7	234-253	25-34
229-248	8	241-260	25-34

	Size	Joint Circum.	Toe Depth	
		\$		
]	1	-	-	
]	2	202-221	21-30	
1	3	208-227	22-31	
4	4	215-234	22-31	
4	5	221-240	23-32	
4	6	228-247	24-33	
4	7	234-253	25-34	
-1	•	0.43 0.00	05.04	

_					
\Box	~	en i	-	~	
	н	ш	150	_	

(Narrower Heel & Shallower Forepart)

	Size	Heel Width	Joint Circum.	Instep Circum.	Toe Depth
D		$\overline{\bullet}$	\$	2	
ā	1	-	-	-	-
2	2	54	214-233	198-217	21-30
E E	3	56	220-239	204-223	22-31
0	4	57	227-246	211-230	22-31
	5	59	233-252	217-236	23-32
	6	60	240-259	224-243	24-33
	7	62	246-265	230-249	25-34
	8	63	253-272	237-256	25-34

	Size	Heel Width	Joint Circum.	Instep Circum.	Toe Depth
		\leftrightarrow	\$	2	
	1	-	-	-	-
ע	2	58	226-245	206-225	21-30
€ 1	3	60	232-251	212-231	22-31
>	4	61	239-258	219-238	22-31
	5	63	245-264	225-244	23-32
	6	64	252-271	232-251	24-33
	7	66	258-277	238-257	25-34
	8	67	265-284	245-264	25-34

	Size	Heel Width	Joint Circum.	Instep Circum.	Toe Depth
		\leftrightarrow	\$	2	
	1	-	-	-	-
	2	50	202-221	190-209	21-30
<u>=</u>	3	52	208-227	196-215	22-31
2	4	53	215-234	203-222	22-31
	5	55	221-240	209-228	23-32
	6	56	228-247	216-235	24-33
	7	58	234-253	222-241	25-34
	8	59	241-260	229-248	25-34

Regular Contour Lasts Features: round toe, 2 x 3mm poron insoles & footbed

	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length		Size	Heel Width	Instep Circ
	⊕	0	\$	2	0	2			\oplus	4
1	55	85	212-224	206-218	24-30	222		1	51	203-21
1	57	87	218-230	212-224	25-31	230		2	53	210-22
1	59	89	224-235	218-230	26-32	238		3	55	216-22
1	61	92	230-242	224-236	26-32	247		4	57	221-23
1	63	94	237-249	231-243	27-33	256		5	59	227-23
1	65	95	243-255	237-249	28-34	265		6	61	232-24
1	67	97	249-261	243-255	29-35	274		7	63	237-24
1	68	100	255-267	250-262	30-36	282		8	64	242-25
1	70	102	260-272	256-268	31-37	291		9	66	248-26
1	71	105	267-279	262-274	31-37	300		10	67	255-26
1	73	107	273-285	268-280	33-39	309		11	69	262-27
1	75	110	280-292	275-287	34-40	318		12	71	268-28
i	76	112	291-303	283-295	35-41	326		13	72	273-28
			25. 000	200 200	00 41	020				

Size	Heel Width	Joint Width	Joint Circum	Instep Circum.	Toe Depth	Last Length
	$\overline{\bullet}$		2			$ \bigcirc $
1	61	91	229-241	219-231	24-30	222
2	63	94	235-247	225-237	25-31	230
3	65	96	241-253	231-243	26-32	238
4	66	99	247-259	237-249	26-32	247
5	68	101	253-265	243-255	27-33	256
6	69	103	259-271	249-261	28-34	265
7	71	105	265-277	255-267	29-35	274
8	72	107	271-283	261-273	30-36	282
9	74	110	277-289	267-279	31-37	291
10	75	113	284-296	274-286	31-37	300
11	77	115	290-302	280-292	33-39	309
12	78	118	296-308	287-299	34-40	318
13	80	120	302-314	293-305	35-41	326

iigui		Width	
\Box		Θ	2
222	1	57	215-227
230	2	59	219-231
238	3	61	225-237
247	4	62	230-242
256	5	64	236-248
265	6	65	242-254
274	7	67	248-260
282	8	68	254-266
291	9	70	260-272
300	10	71	266-278
309	11	73	273-285
318	12	74	279-291
326	13	76	286-298
720			
ast		Heel	

Donna Contour

(Narrower Heel)

el ith	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length		Size	Heel Width	Instep Circun
9	Θ	\$	4					Θ	2
3	79	204-216	199-211	24-30	222		1	49	193-205
4	82	210-222	205-217	25-31	230		2	50	199-211
6	85	216-228	211-223	26-32	238	[3	52	206-218
7	87	222-234	217-229	26-32	247	[4	53	212-224
9	89	228-240	223-235	27-33	256		5	55	218-230
0	92	234-246	229-241	28-34	265		6	56	224-236
2	94	240-252	235-247	29-35	274		7	58	230-242
3	97	246-258	241-253	30-36	282		8	59	236-248
5	99	252-264	247-259	31-37	291		9	61	242-254
6	102	258-270	253-265	31-37	300	[10	62	248-260
В	105	264-276	259-271	33-39	309	[11	64	254-266
0	107	270-282	265-276	34-40	318		12	66	260-272
1	109	276-288	271-283	35-41	326		13	67	266-279

PROFILE LAST FITTING CHARTS

SLIM PROFILE LAST

SIZE	JOINT WIDTH	JOINT CIRC	INSTEP CIRC	TOE DEPTH
1	81	201 - 211	194 - 204	24 - 30
2	83	206 - 216	199 - 209	25 - 31
3	85	211 - 221	204 - 214	26 - 32
4	87	216 - 226	209 - 219	27 - 33
5	89	221 - 231	214 - 224	28 - 34
6	91	226 - 236	219 - 229	29 - 35
7	93	231 - 241	224 - 234	30 - 36
8	95	236 - 246	229 - 239	31 - 37
9	97	241 - 251	234 - 244	32 - 38

STD PROFILE LAST

SIZE	JOINT WIDTH	JOINT CIRC	INSTEP CIRC	TOE DEPTH
1	86	215 - 225	200 - 210	24 - 30
2	88	220 - 230	205 - 215	25 - 31
3	90	225 - 235	210 - 220	26 - 32
4	92	230 - 240	215 - 225	27 - 33
5	94	235 - 245	220 - 230	28 - 34
6	96	240 - 250	225 - 235	29 - 35
7	98	245 - 255	230 - 240	30 - 36
8	100	250 - 260	235 - 245	31 - 37
9	102	255 - 265	240 - 250	32 - 38

WIDE PROFILE LAST

SIZE	JOINT WIDTH	JOINT CIRC	INSTEP CIRC	TOE DEPTH
1	91	222 - 232	203 - 213	24 - 30
2	93	227 - 237	208 - 218	25 - 31
3	95	232 - 242	213 - 223	26 - 32
4	97	237 - 247	218 - 228	27 - 33
5	99	242 - 252	223 - 233	28 - 34
6	101	247 - 257	228 - 238	29 - 35
7	103	252 - 262	233 - 243	30 - 36
8	105	257 - 267	238 - 248	31 - 37
9	105	261 - 271	243 - 253	32 - 38

	Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
D		\leftrightarrow		\$	4		0
	4	62	91	231-249	221-239	31-40	246
ard	5	64	93	237-255	227-245	31-40	254
Stand	6	65	96	244-262	234-252	32-41	263
西	7	67	98	250-268	240-258	33-42	271
S	8	68	101	257-275	247-265	33-42	279
	9	70	103	263-281	253-271	34-43	288
	10	71	106	270-288	260-278	35-44	296
	11	73	108	276-294	266-284	35-44	305
	12	74	111	283-301	273-291	36-45	314
	13	76	113	289-307	279-297	37-46	322
	14	77	116	296-314	286-304	37-46	331

14	77	116	296-314	286-304	37-46	331			
Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length			
	\leftarrow	\oplus	4						
4	66	97	243-261	229-247	31-40	246			
5	68	99	249-267	235-252	31-40	254			
6	69	102	256-274	242-260	32-41	263			
7	71	104	262-280	248-266	33-42	271			
8	72	107	269-287	255-273	33-42	279			
9	74	109	275-293	261-279	34-43	288			
10	75	112	282-300	268-286	35-44	296			
11	77	114	288-306	274-292	35-44	305			
12	78	117	295-313	281-299	36-45	314			
13	80	119	301-319	287-305	37-46	322			
14	81	122	308-326	294-312	37-46	331			

	Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
		\leftrightarrow	θ	\$	4		
	4	58	85	219-237	213-231	31-40	246
	5	60	87	225-243	219-237	31-40	254
	6	61	90	232-250	226-244	32-41	263
	7	63	92	238-256	232-250	33-42	271
	8	64	95	245-263	239-257	33-42	279
	9	66	97	251-269	245-263	34-43	288
	10	67	100	258-276	252-270	35-44	296
1	11	69	102	264-282	258-276	35-44	305
	12	70	105	271-289	265-283	36-45	314
	13	72	107	277-295	271-289	37-46	322
	14	73	110	284-302	278-296	37-46	331

	Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
		\leftrightarrow	θ	4	4		
	4	62	91	243-261	221-239	37-46	246
	5	64	93	249-267	227-245	37-46	254
eeb	6	65	96	256-274	234-252	38-47	263
ě	7	67	98	262-280	240-258	39-48	271
	8	68	101	269-287	247-265	39-48	279
	9	70	103	275-293	253-271	40-49	288
	10	71	106	282-300	260-278	41-50	296
	11	73	108	288-306	266-284	41-50	305
	12	74	111	295-313	273-291	42-51	314
	13	76	113	301-319	279-297	43-52	322
	14	77	116	308-326	286-304	43-52	331

	Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
		\leftrightarrow	Θ	\$	2		
1	4	66	97	255-273	229-247	37-46	246
	5	68	99	261-279	235-253	37-46	254
	6	69	102	268-286	242-260	38-47	263
,	7	71	104	274-292	248-266	39-48	271
	8	72	107	281-299	255-273	39-48	279
	9	74	109	287-305	261-279	40-49	288
	10	75	112	294-312	268-286	41-50	296
	11	77	114	300-318	274-292	41-50	305
	12	78	117	307-325	281-299	42-51	314
	13	80	119	313-331	287-305	43-52	322
	14	81	122	320-338	294-312	43-52	331

Albert

(Narrower Heel)

Size	Heel Width	Instep Circum.
	Θ	2
4	58	211-230
5	60	217-236
6	61	224-243
7	63	230-249
8	64	237-256
9	66	243-262
10	67	250-269
11	69	256-275
12	70	263-282

George (Shallower Forepart)

221-240 227-246 234-253 240-259 247-266 253-272 260-279 266-285 273-292

233-252 239-258

239-258 27-36 246-265 28-37 252-271 29-38 259-278 29-38 265-284 30-39 272-291 31-40 278-297 31-40 285-304 32-41

Size	Heel Width	Instep Circum.
	Θ	2
4	62	219-238
5	64	225-244
6	65	232-251
7	67	238-257
8	68	245-264
9	70	251-270
10	71	258-277
11	73	264-283
12	74	271-290

	Size	Heel Width	IInstep Circum.
		Θ	2
	4	54	203-222
	5	56	209-228
	6	57	216-235
	7	59	222-241
1	8	60	229-248
1 1	9	62	235-254
	10	63	242-261
	11	65	248-267
1 1	12	66	255-274

Size	Heel Width	Ilnstep Circum.	Size	Joint Circum.	Toe Depth
	\oplus	2		A	
4	54	203-222	4	209-228	27-36
5	56	209-228	5	215-234	27-36
6	57	216-235	6	222-241	28-37
7	59	222-241	7	228-247	29-38
8	60	229-248	8	235-254	29-38
9	62	235-254	9	241-260	30-39
10	63	242-261	10	248-267	31-40
11	65	248-267	11	254-273	31-40
12	66	255-274	12	261-280	32-41

Richard (Narrower Heel & Shallower Forepart)

	Size	Heel Width	Joint Circum.	Instep Circum.	Toe Depth
		$\overline{\bullet}$	\$	2	~
o.	4	58	221-240	211-230	27-36
Standard	5	60	227-246	217-236	27-36
2	6	61	234-253	224-243	28-37
to to	7	63	240-259	230-249	29-38
S	8	64	247-266	237-256	29-38
	9	66	253-272	243-262	30-39
	10	67	260-279	250-269	31-40
	11	69	266-285	256-275	31-40
	12	70	273-292	263-282	32-41

	Size	Heel Width	Joint Circum.	Instep Circum.	Toe Depth	
		\leftrightarrow	4	2		
	4	62	233-252	219-236	27-36	
Wide	5	64	239-258	225-244	27-36	
Š	6	65	246-265	232-251	28-37	
_	7	67	252-271	238-257	29-38	
	8	68	259-278	245-264	29-38	
	9	70	265-284	251-270	30-39	
	10	71	272-291	258-277	31-40	
	11	73	278-297	264-283	31-40	
	12	74	285-304	271-290	32-41	
	10 11	71 73	272-291 278-297	258-277 264-283	31-40 31-40	

	Size	Heel Width	Joint Circum.	Instep Circum.	Toe Depth
Slim		$\overline{\bullet}$	A	2	
	4	54	209-228	203-222	27-36
	5	56	215-234	209-228	27-36
	6	57	222-241	216-235	28-37
	7	59	228-247	222-241	29-38
	8	60	235-254	229-248	29-38
	9	62	241-260	235-254	30-39
	10	63	248-267	242-261	31-40
	11	65	254-273	248-267	31-40
	12	66	261-280	255-274	32-41

Regular Contour Lasts Features: round toe, 2 x 3mm poron insoles & footbed

Size	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length
	Θ	0	\$ C	2		
1	55	85	212-225	206-218	24-30	222
2	57	87	218-230	212-224	25-31	230
3	59	89	224-235	218-230	26-32	238
4	61	92	230-242	224-236	26-32	247
5	63	94	237-249	231-243	27-33	256
6	65	95	243-255	237-249	28-34	265
7	67	97	249-261	243-255	29-35	274
8	68	100	255-267	250-262	30-36	282
9	70	102	260-272	256-268	31-37	291
10	71	105	267-279	262-274	31-37	300
11	73	107	273-285	268-280	33-39	309
12	75	110	280-292	275-287	34-40	318
13	76	112	291-303	283-295	35-41	326

	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length	Siz	ze	Heel Width	Instep Circum.
	θ	\$	2		0			$\overline{\bullet}$	2
	91	229-241	213-219	24-30	222	1		57	215-227
1	94	235-247	225-237	25-31	230	_ 2		59	219-231
	96	241-253	231-243	26-32	238	3		61	225-237
1	99	247-259	237-249	26-32	247	4	_	62	230-242
	101	253-265	243-255	27-33	256	_ 5		64	236-248
1	103	259-271	249-261	28-34	265	_ 6)	65	242-254
1	105	265-277	255-267	29-35	274	7		67	248-260
	107	271-283	261-273	30-36	282	- 8	$\overline{}$	68	254-266
	110	277-289	267-279	31-37	291	_ 9		70	260-272
	113	284-296	274-286	31-37	300	10		71	266-278
1	115	290-302	280-292	33-39	309	- 11		73	273-285
	118	296-308	287-299	34-40	318	10		74	279-291
	120	302-314	293-305	35-41	326	- 13	3	76	286-298
_	120	302-314	293-303	30-41	320				

e	Heel Width	Joint Width	Joint Circum.	Instep Circum.	Toe Depth	Last Length		Size	W
	\leftrightarrow	0	4	2		\bigcirc			6
	53	79	204-216	199-211	24-30	222			- 4
	54	82	210-222	205-217	25-31	230		2	
	56	85	216-228	211-223	26-32	238		3	
	57	87	222-234	217-229	26-32	247		4	
	59	89	228-240	223-235	27-33	256		5	
	60	92	234-246	229-241	28-34	265		6	
	62	94	240-252	235-247	29-35	274		7	
	63	97	246-258	241-253	30-36	282		8	
	65	99	252-264	247-259	31-37	291		9	
	66	102	258-270	253-265	31-37	300		10	- (
	68	105	264-276	259-271	33-39	309		11	
	70	107	270-282	265-276	34-40	318		12	
	71	109	276-288	271-283	35-41	326	l l	13	

Αl	bert	Cont	our
(Nar	rrower H	eel)	

Size	Heel Width	Instep Circum.
	0	4
1	51	203-215
2	53	210-222
3 4	55	216-228
4	57	221-233
5	59	227-239
6	61	232-244
7	63	237-249
8	64	242-254
9	66	248-260
10	67	255-267
11	69	262-274
12	71	268-280
13	72	273-285
10	12	213-203

Size	Heel Width	Instep Circum.
	\oplus	2
1	57	215-227
3	59	219-231
3	61	225-237
4	62	230-242
5	64	236-248
6	65	242-254
7	67	248-260
8	68	254-266
9	70	260-272
10	71	266-278
11	73	273-285
12	74	279-291
13	76	286-298

Size	Heel Width	Instep Circum.
	•	2
	49	193-205
2	50	199-211
3	52	206-218
4	53	212-224
5	55	218-230
6	56	224-236
7	58	230-242
8	59	236-248
9	61	242-254
10	62	248-260
11	64	254-266
12	66	260-272
13	67	266-279



COLOURS 4



Modular Footwear Conditions of Sale



Available at no extra cost

- · Half Sizes
- Touch & Close Fastening
- Trial Fitting Service

Buchanan Orthotics

Despatch

Modular footwear will be despatched within seven to ten working days. Shoes requiring adaptions will be despatched within ten to fifteen working days.

As an orthotic company we know the importance of meeting our delivery dates and are always striving to improve on them. On occasion, however, circumstances beyond our control may result in a slight delay and our customer support team will contact you with a revised despatch date should a delay occur.

Carriage

Carriage free to UK mainland.

Payment

Payment is due thirty days from invoice date if account facilities are established. Volume discounts are available for orders of 20, 40 and 80 pairs per month.

20 pairs per month - 5% discount 40 pairs per month - 10% discount 80 pairs per month - 15% discount

Returns

All shoe styles on regular lasts without adaptions/modifications are returnable for credit or exchange if returned within two months in original packaging and in full resaleable condition.

Boot styles are subject to a 30% handling charge.

Split sizes or fittings are non returnable.

Any footwear requested with a through sole will not be returnable unless the through sole is standard for that type of footwear.

Ownership of Goods

The footwear shall remain the absolute property of Buchanan Orthotics Limited, until such time as the customer shall have paid the agreed price in full.

Colours

Our colour swatches shown here may vary from the actual shades due to the limitations of the reproduction processes.