

# LucRo SHOE COLLECTION

Specialist shoes for people  
with diabetes or rheumatism

**LucRo**  
by schein

[www.lucro.de](http://www.lucro.de)



NEW!

LucRo sportic

schein orthopädie  
service

— Germany since 1879 —






Movement improves well-being,  
prevents illness and  
gives us a sense of freedom.

We need to take care of our feet  
if we want to stay happy and free.



— Germany since 1879 —



The logo for LucRo by schein, featuring the brand name in a stylized, outlined font with a red-to-orange gradient. The text "by schein" is written in a smaller, black, sans-serif font below it.

*LucRo*  
by schein

[www.lucro.de](http://www.lucro.de)

*LucRo*  
**SHOE COLLECTION**

[www.lucro.de](http://www.lucro.de)

*LucRo*  
**classic**

*LucRo*  
**kinetic**

*LucRo*  
**sportic**

*LucRo*  
**ergonic**

*Specialist shoes for  
people with diabetes or  
rheumatism*

# *Your feet play a supporting role in your life: give them the protection and attention they deserve.*



Movement is a basic need and is an elementary prerequisite for our health. It increases well-being, prevents diseases and enables us to actively enjoy life. Our feet therefore play a major supporting role for our entire musculoskeletal system.

We usually only notice how important our feet are for our mobility, and thus our well-being, when problems arise – such as pressure sores or minor wounds, which can become a serious health risk for people with diabetes in particular. Rheumatic inflammatory diseases in the foot or toe area are often associated with calluses, misalignments and even toe deformities. These symptoms also require special attention.

With LucRo products, we want to give feet the protection and treatment they deserve. Now with four different shoe collections and a total of 204 items, we offer people with diabetes or rheumatism a variety of options to find their perfect shoe – for their individual foot shape, for their individual symptoms and for their individual taste.

On the street, out in nature or at home, LucRo specialist shoes give you a key piece of your quality of life back: the joy of movement with every step.

**Only the best for your feet. LucRo.**

A handwritten signature in black ink, appearing to read 'A. Schein', located below the bold text.

Arnd Schein  
Schein Orthopädie Service KG

# LucRo SHOE COLLECTION



→ **Indications**

Foot complaints	10–11
Diabetes diagnosis	12–13
The diabetic foot	14–15
Effective foot protection	16–17
Rheumatism diagnosis	18–19
The rheumatic foot	20–21

→ **LucRo specialist shoes**

The foot experts	22–23
New: LucRo sportic	24–27
Collections overview	28–29
LucRo design features	30–33
Our quality promise	34–35
Appendices and care tips	36–37
LucRo accessories	174–179



→ Women's shoes

<b>Sandals</b>	42–46
<b>Shoes with Velcro</b>	47–50
<b>Shoes with laces</b>	51
<b>Boots</b>	52–53
<b>Slippers</b>	54–55

→ Men's shoes

<b>Sandals</b>	60–63
<b>Shoes with Velcro</b>	64–69
<b>Shoes with laces</b>	70–72
<b>Boots</b>	73–75
<b>Slippers</b>	76–77

LucRo classic

→ Women's shoes

<b>Sandals</b>	82–84
<b>Shoes with Velcro</b>	85–88
<b>Shoes with laces</b>	89
<b>Boots</b>	90–91
<b>Slippers</b>	92–93

→ Men's shoes

<b>Sandals</b>	98–100
<b>Shoes with Velcro</b>	101–103
<b>Shoes with laces</b>	104–105
<b>Boots</b>	106–107
<b>Slippers</b>	108–109

LucRo kinetic

→ Women's shoes

<b>Shoes</b>	
<i>with laces/Velcro</i>	114–123
<b>Boots</b>	124–125

→ Men's shoes

<b>Shoes</b>	
<i>with laces/Velcro</i>	130–139
<b>Boots</b>	140–141

LucRo sportic

→ Women's shoes

<b>Sandals</b>	146–148
<b>Shoes with Velcro</b>	149–153
<b>Shoes with stretch leather</b>	152–153
<b>Shoes with laces</b>	153
<b>Boots</b>	154–155
<b>Slippers</b>	156–157

→ Men's shoes

<b>Sandals</b>	162–163
<b>Shoes with Velcro</b>	164–167
<b>Shoes with stretch leather</b>	167
<b>Shoes with laces</b>	168–169
<b>Boots</b>	170–171
<b>Slippers</b>	172–173

LucRo ergonomic

# *Foot complaints* **COMMONPLACE WITH DIABETES OR RHEUMATISM**

## **Millions of people are affected**

Foot complaints are not generally given much attention. With diseases such as diabetes or rheumatism, however, they are anything but rare – millions of people in Germany are affected by diabetic or rheumatic foot complaints. These range from mild calluses through toe deformities up to the dreaded diabetic foot syndrome.



**MORE THAN SIX  
MILLION PEOPLE**

live with a diabetes diagnosis  
in Germany.<sup>1</sup>



**1.5 MILLION  
PEOPLE**

in Germany suffer from  
inflammatory rheumatic illnesses.<sup>2</sup>



**WOMEN ARE  
THREE TIMES  
MORE**

affected than men.<sup>2</sup>



**APPROX. 250,000 PEOPLE SUFFERING FROM DIABETES**

already have open wounds on their feet.<sup>1</sup>



**APPROX. ONE MILLION PEOPLE**

with diabetes are at an increased risk of open foot wounds.<sup>1</sup>



**EVERY FOURTH PERSON WITH DIABETES**

will suffer from diabetic foot syndrome in their life.<sup>1</sup>



**AN ESTIMATED EIGHT TO TEN MILLION PEOPLE**

in Germany suffer from hallux valgus (painful deformity of the big toe).<sup>4</sup>



**90% OF PEOPLE**

with inflammatory rheumatic illnesses will suffer from foot and ankle problems in the course of the disease.<sup>3</sup>

**DIABETES**

PAGES 12–17

**RHEUMATISM**

PAGES 18–21

<sup>1-4</sup> sources: see page 182

# *Diabetes diagnosis* **WALK THROUGH LIFE WITH VIGILANCE**



## **Feet need special attention**

In people with diabetes, pressure sores or blisters caused by shoes that are too narrow or ill-fitting can quickly become a serious health risk. Especially when wounds are noticed much too late or not at all.

Untreated diabetes makes feet more vulnerable. On the one hand, this is due to the impaired metabolism, which makes the skin on the feet more prone to minor cracks and lesions. On the other hand, diabetes is associated with nervous irritation and damage in the long-term, which reduces the sensitivity of the nerves. Injuries are less easy to perceive, if at all. Regular, careful foot checks are therefore essential.

Protective footwear can help to prevent pressure sores in the first place, counteracting inflammatory complications and even amputations at an early stage.

**DIABETES**

*Avoid new wounds*

*Prevent pressure on  
the foot*

*Detect foot lesions  
at an early stage*

**THE RIGHT FOOTWEAR  
CAN HELP**

*Avoid infections and  
amputations*

*Keep moving*

# *The diabetic foot* **KNOW, IDENTIFY AND PREVENT RISKS**



## **What exactly is diabetic foot syndrome?**

Diabetic foot syndrome is one of the most feared secondary diseases of diabetes. Either there is nerve damage in the foot, leading to painless pressure ulcers, or damage to the blood vessels in the foot, causing circulatory disorders. Many cases are present with both forms.

The risk of disease increases with age. Every fourth person with diabetes will suffer from diabetic foot syndrome in the course of their life.

## **How do inflamed foot wounds occur?**

Long-term increased blood sugar levels can damage nerves and blood vessels. The potential consequences are circulatory and sensory disorders.

The triggers for foot infections are usually minor foot injuries or pressure sores, which are not noticed at first due to nerve damage. Pathogens enter the wound and a minor injury becomes a large, inflamed wound. Wound healing is also impaired due to the blood circulation disorder in diabetic feet. Persistent or recurring wounds may be the result. Treatment becomes increasingly difficult. Amputations are required in the worst cases.

## WHAT CAN I DO AS A PATIENT TO PREVENT THESE DREADED CONSEQUENCES?

Prevention is the best protection. People with diabetes can make a decisive contribution to noticing and treating changes to their feet as early as possible:

→ **Daily foot checks:**

Everyone with diabetes should check their feet for changes, preferably daily, so as not to miss any injuries or pressure sores. Family members can of course also help out.

→ **Medical foot clinic:**

In addition to visual checks, regular visits to a medical foot clinic are strongly recommended. The more nerve damage (polyneuropathy) there is, or the more pronounced the symptoms of diabetic foot syndrome are, the more often you should be seen.

→ **Wear the right shoes:**

Well-fitting and, in particular, non-restrictive shoes are more than important as a preventative measure. Orthopaedic shoes for people with diabetes protect the feet with their special features (see also pages 30–33). Clinical studies confirm that the risk of recurring foot ulcers can be reduced by 80% by wearing 'diabetes-appropriate footwear'. It is also advisable to wear suitable footwear with socks at home and not to walk barefoot where possible.



# Effective foot protection

## THE CLINICAL EVIDENCE




### The LucRo classic study

The results of our clinical study carried out at Heinrich Heine University in Düsseldorf, Germany, show that suitable diabetic footwear can make a significant contribution in preventing foot ulcers and/or recurring inflammations (recurrences). It monitored 92 high-risk diabetic patients with healed foot ulcers for recurrence of foot ulcers over a period of 42 months. Some of the patients wore normal footwear, the others wore the LucRo classic specialist diabetes shoe.

### Latest study results on the LucRo ergonomic

A recent study conducted by Dr Dirk Theodor Schraeder, Klinik für Technische Orthopädie, Physikalische Medizin und Rehabilitation (clinic for technical orthopaedics, physical medicine and rehabilitation), at the Hospital zum Heiligen Geist in Geseke, Germany, confirms that the LucRo ergonomic also meets the requirements for a ready-made protective shoe for people with diabetes and is thus a therapeutic tool.





→ Diabetes



# *Rheumatism diagnosis* **DON'T LET IT STOP YOU**



## **Break free of the vicious circle**

Movement improves quality of life and well-being. But if you suffer from rheumatism, you quite often feel the exact opposite: pulling, tearing pain with every movement. Deformations and stiffness in the joints can also aggravate the situation.

Reason enough for many affected persons to reduce their physical activity to the minimum. Yet it's precisely this that increases painful symptoms. It's a vicious circle you must break – with pain relief and the promotion of pain-free movement.

Choosing the right footwear can help a great deal.

**RHEUMATISM**

*Relieve ankles*

*Reduce pain when  
moving*

**THE RIGHT FOOTWEAR  
CAN HELP**

*Bring back the joy  
of movement*

*Provide optimal care for  
foot deformities*

# The rheumatic foot

## REDUCE PAIN, RELIEVE ANKLES




### What is rheumatism exactly?

Rheumatism is a collective concept that includes around 100 different diseases caused by inflammation. The medically correct term is 'rheumatic diseases'. The term 'rheumatism' is Greek in origin and means 'pulling, tearing pain'. This is also precisely what the different forms of the disease have in common: pain and restricted movements.

The most common rheumatic disease is wear-related, degenerative joint disease – primarily arthrosis. There are also the inflammatory joint diseases. Rheumatoid arthritis occurs most frequently. With this disease, the body's defences are misdirected: the immune system fights its own joints and various tissues. Women are affected two or three times more often than men.

### What happens when the feet are affected?

With a rheumatoid disease, the feet and/or the hips and toe joints can also be affected. The most common symptoms of rheumatoid arthritis in the foot are pain, swelling and stiffness, which usually occur bilaterally, i.e. on both ankles. This can impair the gait and foster the development of corns, calluses and toe deformations. Hallux valgus (big toe bunion), hammer toes and claw toes are among the most common defects.





## WHAT CAN I DO AS AN AFFECTED PERSON TO REDUCE PAIN AND RELIEVE THE JOINTS?

Rheumatism has no cure, but there is still a lot we can do to enable people with rheumatism to lead a life as free from symptoms as possible. In addition to medical treatment and surgical interventions, the patients themselves can also help to reduce pain and promote quality of movement.

→ **The right diet:**

A healthy diet agreed with your doctor is an important addition to rheumatoid medications and can help to relieve discomfort. It can also help against obesity, which puts additional strain on joints.

→ **Joint-friendly sports:**

Physical activity has a beneficial effect on the joints of people with rheumatism. It strengthens the muscles, stabilises the gait and soothes pain. Joint-friendly endurance sports such as swimming, cycling, t'ai chi or water aerobics can help.

→ **The right footwear:**

Shoes also have a key part to play in the treatment of rheumatism. Extra-soft materials are used to make orthopaedic specialist shoes, like the LucRo shoe, for people with rheumatism. In combination with their expansive volume, they give toes sufficient freedom of movement. Shock-absorbing soles, support inserts and flattened heels ensure optimal pressure relief and contribute to an improved gait.



# *The foot experts* **FOR PEOPLE WITH DIABETES OR RHEUMATISM**

## **Finding the right footwear together**

Whether you have diabetic foot syndrome or rheumatic discomfort, you don't have to choose the right footwear alone. There are various experts available to help people with diabetes or rheumatism:

### → **The doctor**

who makes the diagnosis checks the health of your feet and prescribes a suitable specialist orthopaedic shoe in the case of an appropriate indication.

### → **The orthopaedic (shoe) technician**

who uses various measurements and in-depth knowledge to select the right shoes for you as an individual, and can also make further necessary adjustments.

### → **The shoe manufacturer**

who makes and develops specialist orthopaedic shoes in close consultation with physicians and orthopaedic shoe technicians, based on current scientific research.



Digital foot scanning system DigiPED

**THE TIP THERM**  
checks the foot's temperature sensitivity.



**Pressure measurements**

In addition to determining the sensitivity of the foot, measuring the pressure distribution of the foot is part of every shoe fitting. Pressure peaks often promote the emergence of wounds under the foot. This makes an even pressure distribution under the sole all the more important.


A simple observation of the soles when standing is done using a podoscope or a digital scanning system. The electronic measurement optimises the quality of individual care: this makes documenting the outcomes of the treatment course easy and sustainable. When making orthopaedic inlays adapted for diabetes, pressure measurement in the shoe is required by your statutory health insurance company.

**Determining foot sensitivity**

**THE TUNING FORK**  
checks depth perception.



**THE MONOFILAMENT**  
checks surface perception.



Pressure under the sole

**Measuring the pressure distribution of the foot for the perfect shoe fit**



*Out into nature,*  
**INTO A FULL LIFE**

The new LucRo collection



→ New LucRo sportic collection

**PERFECTLY  
EQUIPPED**

for strolls and walks  
as well as other  
leisure activities



*Stability on slippery  
surfaces*

*Diabetes-friendly  
foot protection*

*LucRo*  
**sportic**

*Sturdy functional and  
leisure shoe*

*Athletic fit*

*Trendy colours*

**NEW!**

*LucRo sportic*

# *LucRo sportic* **ACTIVE EVERY DAY**



## **Movement is the best medicine**

Movement is not only healthy, it is also beneficial. In people with diabetes, it helps to reduce blood sugar levels. People with rheumatism can stop or slow down the progress of joint wear through regular movement.

But just as important: physical activity also promotes emotional well-being, and therefore overall quality of life. Being physically active doesn't have to be in a sports club or gym. It can be a long walk, a group trip into the countryside or a forest hike. The LucRo sportic is perfect for just this.



## LUCRO SPORTIC: FUNCTIONAL AND SPECIALIST SHOE IN ONE

Just like the other three LucRo shoe series classic, kinetic and ergonic, the latest collection also ensures perfect all-round protection that takes the pressure off the foot and reduces the risk of recurrent diabetic foot syndrome.

→ **Anti-slip and weatherproof:**

With its sturdy sole, its weather-resistant material and athletic fit, the new LucRo sportic is the ideal outdoor shoe alternative to your solid everyday shoe. On gravel paths, hiking trails or cross-country – it gives feet a secure hold without pressure.

→ **Innovative 'rebound effect':**

A special material in the midsole – an innovative particle foam, which creates what is known as the rebound effect in the latest running shoe models from well-known sporting-goods manufacturers, is also used in the new LucRo sportic. The special feature of the new sole material: if a high level of pressure is created when stepping, the sole immediately springs back to its original form. At low pressure, the material cushions the foot and absorbs the energy. This means the shoe wearer needs much less energy when walking and the joints are relieved.

*Collection pages 110–141*

# LucRo

## THE BEST FOR EVERY FOOT



### The classic one

**The LucRo classic can protect against recurring ulcers in diabetic foot syndrome with almost 80% certainty.**

This was clinically proven in a 42-month study (see also page 16). But pressure relief and protection against recurring ulcers are certainly not everything the LucRo classic has to offer.

#### Perfect all-round protection plus:

- anatomically shaped, classically formed last
- large and diverse selection of 72 models

*Pages 38–77*



### The slimmer one

**The LucRo kinetic is the slimmer-cut alternative that also provides ulcer-recurrence protection.**

The feature of the LucRo kinetic line is a modified fit. We have given the last form a slimmer and narrower rear foot in the heel. This gives the shoe a different look and ensures an optimised heel guide especially for slimmer heels.

#### Perfect all-round protection plus:

- last in a narrower natural form
- special kineticBags in the soles that relieve the joints
- large and diverse selection of 46 models

*Pages 78–109*





## The sporty one

**New to the range: the specialist shoe for people with diabetes who want more active leisure time.**

The latest LucRo collection also has the protective properties for diabetes-friendly foot care. Thanks to its anti-slip sole and sturdy design, it ensures stability especially over rough terrain.

### Perfect all-round protection plus:

- last in a narrower natural form
- anti-slip, sturdy TPU outsole
- innovative, joint-protecting midsole
- large and diverse selection of 48 models

Pages 110–141



## The adaptable one

**In addition to proven ulcer-recurrence protection, the LucRo ergonomic also provides ideal features for rheumatism care.**

A recent diabetes study shows that the LucRo ergonomic effectively protects against recurring foot ulcers. Its joint-friendly features also reduce the symptoms of rheumatism patients. This is why the LucRo specialist shoe is listed for both indications.

### Perfect all-round protection plus:

- last in a natural form
- stiffening sole for optimum pressure distribution (must be firmly glued in for diabetes!)
- large and diverse selection of 38 models

Pages 142–173



# LucRo features

## GOOD DESIGN ON THE OUTSIDE, GOOD PROTECTION ON THE INSIDE



### Trendy designs just like a modern everyday shoe

Classic or athletic fit, narrow foot shape or expansive volume – in addition to optimum all-round protection, all LucRo collections have another thing in common: their appealing design. The therapeutic benefit that comes with each LucRo shoe is packaged in a modern-looking shoe perfect for every day.

### Interchangeable stiffening soles for LucRo ergonomic


Depending on the indication and individual form of your complaints, you can choose from up to three different stiffening soles. Note: In the case of diabetes, the stiffening sole must be glued in firmly!



Three interchangeable  
stiffening soles

### Breathable Sympatex membrane for weather-resistant LucRo boots

The use of a breathable Sympatex membrane in LucRo boots makes them 100% waterproof. This intelligent combination of materials guarantees freedom of movement in all weather conditions. It is made from recyclable, harmless polyether-ester and is therefore absolutely environmentally and skin-friendly.



100%  
waterproof





### Joint-friendly kineticBags in the LucRo kinetic

The kineticBags specially incorporated in the outsole absorb the impact energy and ensure additional joint relief.

### Midsole with 'rebound effect' in the LucRo sportic

A special midsole material made of expanded thermoplastic polyurethane (TPU) ensures comfort and excellent running characteristics. The sole immediately rebounds into its original shape after stepping. Running therefore requires less energy and joints are relieved: the so-called rebound effect.



### Soft, high-quality upper materials

We use sophisticated and high-quality materials from sustainable production (see also pages 34–35). Selected LucRo ergonomic models are also available with stretch leather, which offers people with rheumatism in particular additional comfort.



### Wide opening makes it easy to put them on

No pressure should be applied to the foot, even when putting on the shoes. Thanks to its wide opening, it's easy, comfortable and safe to step into the LucRo shoe.



# *LucRo features* **PERFECTION DOWN TO THE LAST DETAIL**







→ LucRo design features

**Four different widths**


The four widths S, M, W and +10 (10 mm additional upper height) offer a good selection of care options in each LucRo shoe collection.



Three-part, shock-absorbing insole made from special soft foam




Different firmness in front and rear foot




Stiffened insole with LucRo classic, kinetic and sportic

Stiffening element



Relieving sole with heel and forefoot roll for controlled rolling



LucRo outer sole for even pressure distribution under the foot

Different outer soles depending on LucRo collection

# *Our quality promise* **SUSTAINABLE AT EVERY STEP**

## **The highest standards in material and production**


As a responsible family business, it is a matter of course for us to produce exemplary quality. Our high quality standards don't just apply to our products, but to the entire business. The issue of sustainability also plays a decisive role. This applies to our own manufacturing processes as well as to the choice of our material suppliers. We only work with qualified manufacturers who meet our high standards of high-quality processing and sustainable production.

## **TÜV-certified quality management**

We request an independent, external quality inspection every year. Since 1997, our company's performance has been distinguished by the TÜV Rheinland quality-management certificate. This certifies that our work processes are systematised and the manufacture of high-quality products is ensured.

## **QUALITY YOU CAN TRUST**

- **artisanal**  
Excellent quality for over 130 years
- **sustainable**  
Sustainability is a key quality feature in our own processes and supplier selection
- **medical and scientific**  
Medical efficacy of LucRo classic and ergonomic demonstrated by clinical studies
- **independently tested**  
Independent and external quality inspection as well as quality management certified by TÜV Rheinland





## ENVIRONMENTALLY FRIENDLY

LucRo leather – beautifully worked, sustainably sourced

### Environmentally friendly leather production

Tannins are an important part of preserving animal skins and making them into leather. These days, chemically produced tanning agents are generally used to meet modern needs. The key here is correct implementation, taking account of all environmental protection requirements. This applies to the tanning process as well as to the disposal of chemical substances.

We only use appropriate and sustainably produced leather types so you can always be sure of the highest quality at LucRo.

This symbol tells you which LucRo shoes are made from pure vegetable-tanned leather.



### The alternative: purely vegetable-based tanning

One alternative to chemical tanning is vegetable tanning. This uses only tannins from nature, such as extracts of wood or rhubarb. No artificial or synthetic additives are used. This process is one of the healthiest for humans and the environment. Vegetable-tanned leather also protects skin flora. This process costs the tanner much more effort and time: the entire manufacturing process up to the finished leather can often take up to two years. The residues from vegetable tannins are easy to recycle and are used as agricultural fertilisers.

# Foot protection at home

## COMFORT THAT GOES FURTHER

### LucRo slipper models

You are always moving, especially at home. People with diabetes should not go barefoot here either, and should ensure adequate foot protection. We bring the LucRo comfort you're used to into the home with our various slipper styles.

Easy to step into, lightweight and made with soft materials – you'll want to wear them and feel good all the time.

### LucRo socks for even more protection

Our special LucRo socks have all the features you need for optimum foot protection. Unlike normal socks, they are padded across the entire sole area. As they are hand-knotted, no pressure points occur in the toe area and the elastic-free soft cuff also does not exert any pressure. The material composition is particularly skin-friendly. Silver fibres in the cotton also have an antibacterial effect.

### LucRo Komfort diabetic socks

with all the features for optimum foot protection



## IMPORTANT LUCRO CARE TIPS

Only a fully functional shoe provides optimum protection. With a little care and careful handling, you can retain that special LucRo quality for a long time.

→ **Give your shoes a rest:**

Changing your shoes after a single wear and giving them a rest for 24 hours significantly extends the life of your LucRo shoes.

→ **Change your insoles regularly:**

They can lose their functionality after a time due to high stress.

→ **Use a shoehorn:**

This means you'll protect the extended rear cap and get the optimum heel guide.

→ **Let your shoes dry at room temperature:**

Slow drying (never place them in front of a heater) and regular care with a high-quality shoe polish or nubuck stick protects the leather from drying out.

→ **Protect your shoes before the first wear:**

By spraying the upper material with a protective spray (preferably regularly), the shoes become better protected from moisture.

→ **Regularly brush the leather:**

Regular care with leather brushes and care creams (LucRo care accessories, see page 176) means you will maintain the high quality of the leather.



**Anatomically shaped**  
classically formed last

**Seamless, padded**  
microfibre lining

**Four different  
widths**

→ LucRo classic Damen

classic WOMEN



LucRo  
**classic**

**FOR  
HER**



*Monika*

Pages 43–44



*Anita*

Page 45



*Marisa*

Page 46



*Sylvia*

Page 48



*Yvonne*

Page 49



*Ilona*

Page 50



*Andrea*

Page 51



*Doris*

Page 53



*Fiona*

Page 55



*Hildegard*

Page 55



*Frieda*

Page 55



*LucRo classic*  
**WOMEN'S COLLECTION**



***Sandals***

From page 42

***Shoes***  
*with Velcro/laces*

From page 47

***Boots***

From page 52

***Slippers***

From page 54



*LucRo classic*  
**SANDALS**



**LucRo classic MONIKA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22320	calfskin nappa/patent	black/bordeaux	● ●	34-45	S, M, W	upper
22321	calfskin nappa	white/beige	● ●	34-45	S, M, W	upper
22325	nubuck	beige/cognac	● ●	34-45	S, M, W	upper



22425

Item no.



22328

Item no.



Monika  
Model

22423

Item no.

### LucRo classic **MONIKA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22328	nubuck	ocean	●	34-45	S, M, W	upper
22423	goatskin nappa metallic	bronze/sand	● ●	36-43	S, M, W	upper
22425	goatskin nappa metallic	sand/bronze	● ●	36-43	S, M, W	upper



Anita  
Model  
22409  
Item no.



22403  
Item no.



22401  
Item no.



22400  
Item no.

LucRo classic ANITA

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22400	calfskin nappa	black	●	36-43	M, W	S, upper
22401	calfskin nappa	pearl	●	36-43	M	S, W, upper
22403	antiqued calfskin nappa	brown	●	36-43	M	S, W, upper
22409	goatskin nappa metallic	anthracite	●	36-43	M	S, W, upper



KIND VEGETABLE  
LEATHER TANNING



22063  
Item no.



22066  
Item no.



22065  
Item no.

### LucRo classic **MARISA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22063	goatskin nappa metallic	bronze		36-43	M	S, W, +10, upper
22065	nubuck	camel		36-43	M, W	S, +10, upper
22066	vegetable nappa	red		36-43	M	S, W, +10, upper



*LucRo classic*  
**SHOES**

→ LucRo classic women's shoes with Velcro



### LucRo classic SYLVIA

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22040	buckskin/nubuck	black	●	34-45	S, M, W, +10, upper	/
22046	vegetable nappa	red	●	36-43	M	S, W, +10, upper
22049	antiqued calfskin nappa	grey	●	36-43	M, W, +10	S, upper
22445	goatskin nappa metallic	sand	●	36-43	M	S, W, +10, upper





22301

Item no.



Yvonne

Model



22300

Item no.

LucRo classic **YVONNE**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22300	calfskin nappa/nubuck/ patent	black/bordeaux	● ●	36-43	M, W	S, +10, upper
22301	calfskin nappa	white/beige	○ ●	36-43	M	S, W, +10, upper

→ LucRo classic women's shoes with Velcro



22362  
Item no.



22363  
Item no.



22360  
Item no.

*Ilona*  
Model

### LucRo classic **ILONA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22360	calfskin nappa/patent	black/bordeaux	● ●	36-43	M, W, +10	S, upper
22362	vegetable nappa	black	●	34-45	M, W, +10	S, upper
22363	goatskin nappa metallic	bronze/sand	● ●	36-43	M	S, W, +10, upper

→ LucRo classic women's shoes with laces

22190

Item no.



22495

Item no.



Andrea

Model

LucRo classic **ANDREA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22190	calfskin nappa/patent	black/bordeaux	● ●	36-43	M, W	S, +10, upper
22495	goatskin nappa metallic	sand	●	36-43	M	S, W, +10, upper



*LucRo classic*  
**BOOTS**



### LucRo classic **DORIS**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22150	calfskin nappa Hydro	black	●	36-43	M, W	S, +10, upper
22159	goatskin nappa metallic/ fur	anthracite/grey	● ●	36-43	M	S, +10, upper



*LucRo classic*  
**SLIPPERS**

→ LucRo classic women's slippers



22399  
Item no.



22394  
Item no.

LucRo classic **HILDEGARD**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22394	felt	beige	●	36-43	S, M, W	upper
22399	felt	grey	●	36-43	S, M, W	upper



22413  
Item no.



22419  
Item no.

LucRo classic **FIONA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22413	felt/ antiqued calfskin nappa	beige	●	36-43	M, W	S, W, upper
22419	felt/nubuck	grey/black	● ●	36-43	M, W	S, W, upper

22389  
Item no.



22384  
Item no.

LucRo classic **FRIEDA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
22384	felt/nubuck	beige	●	36-43	M, W	S, upper
22389	felt/nubuck	grey	●	36-43	M, W	S, upper



**Anatomically shaped**  
classically formed last

**Seamless, padded**  
microfibre lining

**Four different  
widths**

**Lucho**



→ LucRo classic men



classic | MEN

LucRo  
**classic**

**FOR  
HIM**



*Eric*

Page 61



*Achim*

Page 62



*Peter*

Page 62



*Oliver*

Page 63



*Paul*

Page 63



*Harald*

Page 65



*Wolfgang*

Page 66



*Volker*

Page 66



*Thorsten*

Page 67



*Walter*

Page 68



*Sven*

Page 69



*Philip*

Page 70



*William*

Page 71



*Björn*

Page 72



*Lars*

Page 72



*Jens*

Page 74



*Frank*

Page 75



*Kurt*

Page 75



*Heinrich*

Page 77



*Fritz*

Page 77



*Falco*

Page 77

*LucRo classic*  
**MEN'S COLLECTION**



***Sandals***

From page 60

***Shoes***  
*with Velcro/laces*

From page 64

***Boots***

From page 73

***Slippers***

From page 76



*LucRo classic*  
**SANDALS**

→ LucRo classic men's sandals



21273  
Item no.



21270  
Item no.



21278  
Item no.



21276  
Item no.

LucRo classic **ERIC**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21270	nubuck	black	●	37-50	S, M, W	upper
21273	nubuck	brasil	●	37-50	S, M, W	upper
21276	calfskin nappa	cognac	●	40-47	S, M, W	upper
21278	vegetable nappa	grey	●	40-47	S, M, W	upper

→ LucRo classic men's sandals



### LucRo classic **ACHIM**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21190	calfskin nappa	black	●	40-47	M, W, upper	S
21192	calfskin nappa	cognac	●	40-47	M, W	S, upper



### LucRo classic **PETER**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21200	calfskin nappa	black	●	40-47	M, W, upper	S, +10



Paul  
Model


21240  
Item no.

LucRo classic **PAUL**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21240	calfskin nappa	black	●	40-47	M, W, +10	S, upper



21090  
Item no.



21097  
Item no.

Oliver  
Model



21096  
Item no.



21095  
Item no.

LucRo classic **OLIVER**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21090	calfskin nappa	black	●	40-47	M, W +10	S, upper
21095	calfskin nappa	cognac	●	40-47	M, W, +10	S, upper
21096	nubuck	brasil	●	40-47	M, W, upper	S, +10
21097	nubuck	camel	●	40-47	M	S, W, +10, upper



*LucRo classic*  
**SHOES**



→ LucRo classic men's shoes with Velcro



21360  
Item no.



21365  
Item no.

### LucRo classic HARALD

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21360	calfskin nappa/nubuck/fabric	black	●	40-47	M, W	S, upper
21365	calfskin nappa/nubuck/fabric	camel/desert/mottled brown	● ● ●	40-47	M, W	S, upper

→ LucRo classic men's shoes with Velcro



### LucRo classic **WOLFGANG**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21440	calfskin nappa	black	●	40-47	M	S, W, upper
21443	calfskin nappa	cognac	●	40-47	M	S, W, upper



### LucRo classic **VOLKER**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21330	buckskin	black	●	37-50	M, W, +10	S, upper
21333	nubuck	brasil	●	40-47	M, W	S, +10, upper

*Thorsten*  
Model

21020  
Item no.



21029  
Item no.



21025  
Item no.



**LucRo classic THORSTEN**


Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21020	calfskin nappa/nubuck	black	●	37-50	S, M, W, +10, upper	/
21025	calfskin nappa	cognac	●	40-47	M, W	S, +10, upper
21029	antiqued calfskin nappa	grey	●	40-47	M, W	S, +10, upper

→ LucRo classic men's shoes with Velcro

21370  
Item no.



KIND VEGETABLE  
LEATHER TANNING



KIND VEGETABLE  
LEATHER TANNING



21373  
Item no.

Walter  
Model



KIND VEGETABLE  
LEATHER TANNING



21379  
Item no.

LucRo classic **WALTER**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21370	vegetable nappa	black	●	37-50	M, W, +10	S, upper
21373	vegetable nappa	brown	●	40-47	M	S, W, +10, upper
21379	vegetable nappa	grey	●	40-47	M	S, W, +10, upper



21050

Item no.

Sven  
Model



21053

Item no.



LucRo classic **SVEN**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21050	buckskin	black	●	37-50	S, M, W, +10, upper	/
21053	buckskin	chestnut	●	37-50	S, M, W, +10, upper	/

→ LucRo classic men's shoes with laces



### LucRo classic PHILIP

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21180	calfskin nappa	black	●	37-50	M, W	S, +10, upper
21185	calfskin nappa	cognac	●	40-47	M	S, W, +10, upper



21280  
Item no.



William  
Model



21285  
Item no.

LucRo classic **WILLIAM**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21280	calfskin nappa	black	●	37-50	M, W, +10, upper	S
21285	calfskin nappa	cognac	●	40-47	M	S, W, +10, upper

→ LucRo classic men's shoes with laces

21030

Item no.



*Björn*  
Model

### LucRo classic **BJÖRN**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21030	calfskin nappa	black	●	37-50	M, W	S, +10, upper

21040

Item no.



*Lars*  
Model

### LucRo classic **LARS**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21040	buckskin	black	●	40-47	M, W, upper	S, +10





classic MEN

*LucRo classic*  
**BOOTS**

→ LucRo classic men's boots with Velcro



### LucRo classic JENS

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21260	calfskin nappa Hydro	black	●	37-50	M, W, +10	S, upper
21263	nubuck Hydro	brasil	●	40-47	M	S, W, +10, upper



21320

Item no.



Frank  
Model

### LucRo classic FRANK

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21320	calfskin nappa Hydro with warm lining	black	●	37-50	M, W, +10	S, upper

→ LucRo classic men's boots with laces



21310

Item no.



KIND VEGETABLE  
LEATHER TANNING



Kurt  
Model

21313

Item no.



### LucRo classic KURT

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
21310	vegetable nappa	black	●	40-47	M	S, W, +10, upper
21313	nubuck Hydro	brasil	●	40-47	M	S, W, +10, upper



*LucRo classic*  
**SLIPPERS**

→ LucRo classic men's slippers

21399  
Item no.



Heinrich  
Model

### LucRo classic HEINRICH

Item no.	Material	Colour	Size	Width ex stock	Delivery time four weeks
21399	felt	grey	40-47	S, M, W	upper

Fritz  
Model

21389  
Item no.



### LucRo classic FRITZ

Item no.	Material	Colour	Size	Width ex stock	Delivery time four weeks
21389	felt/nubuck	grey	40-47	S, M, W	upper

Falko  
Model

21409  
Item no.



### LucRo classic FALKO

Item no.	Material	Colour	Size	Width ex stock	Delivery time four weeks
21409	felt/nubuck	grey/black	40-47	M, W	S, upper



**Anatomically shaped**  
narrow, natural-form last

**Seamless, padded**  
microfibre lining

**Special kineticBags**  
for additional joint  
relief

→ LucRo kinetic women



LucRo  
**kinetic**

**FOR  
HER**

kinetic WOMEN



*Ingrid*

Page 83



*Kerstin*

Page 84



*Brigitte*

Page 86



*Sabine*

Page 87



*Elke*

Page 88



*Gabriele*

Page 89



*Eva*

Page 89



*Angelika*

Page 91



*Franziska*

Page 93



*LucRo kinetic*  
**WOMEN'S COLLECTION**



***Sandals***

From page 82

***Shoes***  
*with Velcro/laces*

From page 85

***Boots***

From page 90

***Slippers***

From page 92



*LucRo kinetic*  
**SANDALS**

→ LucRo kinetic women's sandals



25055

Item no.



KIND VEGETABLE  
LEATHER TANNING



25056

Item no.



25059

Item no.



25050

Item no.



*Ingrid*  
Model

LucRo kinetic **INGRID**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25050	calfskin nappa	black	●	34-45	S, M, W	upper
25055	nubuck	desert	●	36-43	S, M, W	upper
25056	vegetable nappa	red	●	36-43	S, M, W	upper
25059	goatskin nappa metallic/ embossed snakeskin	anthracite/grey	● ●	36-43	S, M, W	upper



**LucRo kinetic KERSTIN**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25350	calfskin nappa	black/anthracite	● ●	36-43	S, M, W	upper
25351	calfskin nappa/patent	pearl/bordeaux	● ●	36-43	S, M, W	upper
25353	goatskin nappa metallic	bronze/sand	● ●	36-43	S, M, W	upper








kinetic WOMEN

*LucRo kinetic*  
**SHOES**

→ LucRo kinetic women's shoes with Velcro



KIND VEGETABLE  
LEATHER TANNING



### LucRo kinetic **BRIGITTE**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25020	calfskin nappa	black	●	36-43	M, W	S, upper
25023	goatskin nappa metallic	bronze	●	36-43	M	S, W, upper
25026	vegetable nappa	red	●	36-43	M	S, W, upper
25029	goatskin nappa metallic	anthracite	●	36-43	M	S, W, upper

25070

Item no.



25379

Item no.



Sabine

Model

25375

Item no.



LucRo kinetic **SABINE**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25070	calfskin nappa/patent	black/bordeaux	● ●	34-45	S, M, W, +10	upper
25375	goatskin nappa metallic	sand/bronze	● ●	36-43	M	S, W, +10, upper
25379	goatskin nappa metallic/ embossed snakeskin	anthracite/grey	● ●	36-43	M	S, W, +10, upper

→ LucRo kinetic women's shoes with Velcro



25112

Item no.



25110

Item no.



25117

Item no.



25313

Item no.

## LucRo kinetic ELKE

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25110	calfskin nappa/nubuck	black	●	36-43	M, +10	S, W, upper
25112	antiqued calfskin nappa	taupe/white	● ○	36-43	M, W	S, +10, upper
25117	vegetable nappa	black	●	36-43	M, W	S, +10, upper
25313	goatskin nappa metallic	bronze	●	36-43	M	S, W, +10, upper



→ LucRo kinetic women's shoes with laces

25030

Item no.



Gabriele  
Model

### LucRo kinetic **GABRIELE**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25030	calfskin nappa	black	●	36-43	S, M, W	+10, upper

Eva

Model

25100

Item no.



25309

Item no.



25102

Item no.



### LucRo kinetic **EVA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25100	antiqued calfskin nappa	black	●	36-43	M	S, W, +10, upper
25102	antiqued calfskin nappa	taupe/white	● ○	36-43	M, W	S, +10, upper
25309	goatskin nappa metallic	anthracite	●	36-43	M	S, W, +10, upper





*LucRo kinetic*  
**BOOTS**

→ LucRo kinetic women's boots with Velcro



 **Hydro**  
water-resistant leather

**LucRo kinetic ANGELIKA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25010	calfskin Hydro/patent	black/bordeaux		36-43	M, W, +10	S, upper
25413	goatskin nappa metallic/ fur	bronze/brown		36-43	M	S, W, +10, upper

kinetic WOMEN



*LucRo kinetic*  
**SLIPPERS**

→ LucRo kinetic women's slippers



**LucRo Komfort diabetic socks**  
for even more cosy comfort at home



**25143**  
Item no.

*Franziska*  
Model

**25149**  
Item no.



**LucRo kinetic FRANZISKA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
25143	felt/ antiqued calfskin nappa	beige	●	36-43	M, W	S, upper
25149	felt/nubuck	grey/black	● ●	36-43	M, W	S, upper

**Anatomically shaped**  
narrow, natural-form last

**Seamless, padded**  
microfibre lining

**Special kineticBags**  
for additional joint  
relief

→ LucRo kinetic men



LucRo  
**kinetic**

**FOR  
HIM**

Sympatex®

kinetic MEN



*Klaus*

Page 99



*Hans*

Page 100



*Rolf*

Page 102



*Gerhard*

Page 103



*Norbert*

Page 104



*Holger*

Page 105



*Rüdiger*

Page 107



*Franz*

Page 109



*LucRo kinetic*  
**MEN'S COLLECTION**



***Sandals***

From page 98

***Shoes***  
*with Velcro/laces*

From page 101

***Boots***

From page 106

***Slippers***

From page 108



*Hans*

Model

24013

Item no.

*LucRo kinetic*  
**SANDALS**



24020

Item no.

*Klaus*  
Model



24323

Item no.



24025

Item no.

## LucRo kinetic **KLAUS**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
24020	calfskin nappa	black	●	40-47	M, W	S, +10, upper
24025	nubuck	camel	●	40-47	M	S, W, +10, upper
24323	calfskin nappa	cognac	●	40-47	M	S, W, +10, upper



### LucRo kinetic HANS

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
24010	nubuck	black/grey	● ●	37-50	S, M, W	upper
24013	calfskin nappa	cognac	●	40-47	S, M, W	upper
24019	antiqued calfskin nappa/ nubuck	grey	●	40-47	S, M, W	upper



kinetic MEN

*LucRo kinetic*  
**SHOES**



→ LucRo kinetic men's shoes with Velcro



24052  
Item no.



24053  
Item no.



24453  
Item no.



24050  
Item no.



### LucRo kinetic **ROLF**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
24050	calfskin nappa	black	●	37-50	S, M, W, +10	upper
24052	vegetable nappa	black	●	37-50	M, W, +10	S, upper
24053	calfskin nappa	dark brown	●	40-47	M, W, +10	S, upper
24453	calfskin nappa	cognac	●	40-47	M	S, W, +10, upper



*Gerhard*  
Model



24099



Item no.



24093

Item no.

*LucRo kinetic* **GERHARD**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
24093	calfskin nappa/nubuck	cognac		40-47	M, +10	S, W, upper
24099	antiqued calfskin nappa	grey		40-47	M, W	S, +10, upper

→ LucRo kinetic men's shoes with laces



24043  
Item no.

24040  
Item no.



24042  
Item no.



KIND VEGETABLE  
LEATHER TANNING

## LucRo kinetic **NORBERT**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
24040	nubuck/fabric	black	●	40-47	M, W	S, +10, upper
24042	vegetable nappa/ mesh	black/grey	● ●	40-47	M	S, W, +10, upper
24043	nubuck/fabric	brasil/mottled brown	● ●	40-47	M, W	S, +10, upper





24030

Item no.



24034

Item no.



*Holger*  
Model

24033

Item no.

## LucRo kinetic **HOLGER**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
24030	calfskin nappa	black	●	37-50	S, M, W, +10	upper
24033	calfskin nappa	dark brown	●	40-47	M	S, W, +10, upper
24034	calfskin nappa	cognac	●	40-47	M	S, W, +10, upper



*LucRo kinetic*  
**BOOTS**



24063

Item no.



24060

Item no.



## LucRo kinetic RÜDIGER

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
24060	calfskin nappa Hydro	black	●	40-47	M, W, +10	S, upper
24063	nubuck Hydro	brasil	●	40-47	M	S, W, +10, upper



*LucRo kinetic*  
**SLIPPERS**

**LucRo Komfort diabetic socks**

for even more cosy comfort at home



**24119**  
Item no.

*Franz*  
Model



*LucRo kinetic* **FRANZ**

Item no.	Material	Colour	Size	Width ex stock	Delivery time four weeks
24119	felt/nubuck	grey/black	40-47	M, W	S, upper



**Innovative midsole**

with joint-friendly  
'rebound effect'

**Anti-slip**

TPU outsole

**Four different  
widths**

**Sporty design**

in trendy colours

→ LucRo sportic women



LucRo  
**sportic**

**FOR  
HER**



*Maren*

Pages 115–117



*Nina*

Pages 118–119



*Bettina*

Page 120



*Anke*

Pages 121–123



*Bianka*

Page 125





*LucRo sportic*  
**WOMEN'S COLLECTION**



**Shoes**  
*with laces/Velcro*

From page 114

**Boots**

From page 124



*LucRo sportic*  
**SHOES**

→ LucRo sportic women's shoes with laces



30020  
Item no.



30021  
Item no.

Includes spare pair of laces:  
LucRo flat laces in white, 110 cm



*Maren*  
Model

30024  
Item no.

Includes spare pair of laces:  
LucRo flat laces in white,  
110 cm

## LucRo sportic MAREN

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30020	calfskin nappa/ patent fabric	black	●	36-43	M, W	S, +10, upper
30021	calfskin nappa	white	○	36-43	M, W	S, +10, upper
30024	microfibre/ patent canvas	dark brown/brown	● ●	36-43	M, W	S, +10, upper

sportic WOMEN



**30023**  
Item no.



**30022**  
Item no.



**30025**  
Item no.



### LucRo sportic **MAREN**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30022	calfskin nappa/ patent canvas	beige	●	36-43	M, W	S, +10, upper
30023	microfibre/ leopard embossed nubuck	black/white	● ○	36-43	M, W	S, +10, upper
30025	goatskin nappa metallic/ patent fabric	sand/silver	● ●	36-43	M	S, W, +10, upper

30026

Item no.



30027

Item no.




30028

Item no.



**Includes spare pair of laces:**  
LucRo flat laces in dark blue,  
110 cm



## LucRo sportic MAREN

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30026	microfibre/ 3D mesh	red/white	<span style="color: red;">●</span> <span style="border: 1px solid black; border-radius: 50%; padding: 2px;">○</span>	36-43	M	S, W, +10, upper
30027	nubuck/calfskin nappa/ 3D mesh	grey/white	<span style="color: grey;">●</span> <span style="border: 1px solid black; border-radius: 50%; padding: 2px;">○</span>	36-43	S, M, W, +10	upper
30028	microfibre/ patent nubuck	dark blue	<span style="color: blue;">●</span>	36-43	S, M, W, +10	upper

→ LucRo sportic women's shoes with laces



*Nina*  
Model

30037  
Item no.



30038  
Item no.



Includes spare pair of laces:  
LucRo flat laces in blue, 130 cm



### LucRo sportic **NINA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30037	microfibre/ 3D mesh	blue/light blue	<span style="color: #0070C0;">●</span> <span style="color: #ADD8E6;">●</span>	36-43	M, W	S, +10, upper
30038	patent nubuck/ patent canvas/microfibre	dark blue	<span style="color: #003366;">●</span>	36-43	M, W	S, +10, upper



30030  
Item no.



30039  
Item no.



30032  
Item no.



Includes spare pair of laces:  
LucRo flat laces in white, 130 cm

LucRo sportic **NINA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30030	microfibre/ 3D mesh	black	●	36-43	S, M, W, +10	upper
30032	patent nubuck	bordeaux	●	36-43	M, W, +10	S, upper
30039	embossed patent/ microfibre	taupe/grey	● ●	36-43	M	S, W, +10, upper



## LucRo sportic **BETTINA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30050	vegetable nappa/ microfibre	black	●	36-43	M, W	S, +10, upper
30051	calfskin nappa/ 3D mesh	white/blue	○ ●	36-43	M	S, W, +10, upper
30052	vegetable nappa/ microfibre	grey/dark grey	● ●	36-43	S, M, W, +10	upper





Anke  
Model

30010  
Item no.

30011  
Item no.



30012  
Item no.

LucRo sportic ANKE

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30010	embossed patent/ calfskin nappa/microfibre	black	●	36-43	S, M, W, +10	upper
30011	calfskin nappa/microfibre	white/pink	○ ●	36-43	M	S, W, +10, upper
30012	calfskin nappa/ leopard embossed nubuck	porcelain	●	36-43	M, W	S, +10, upper



## LucRo sportic ANKE

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30015	goatskin nappa metallic/ patent fabric	sand/silver		36-43	M, W	S, +10, upper
30017	calfskin nappa/3D mesh/ microfibre	white/grey/petrol		36-43	S, M, W, +10	upper



**30018**  
Item no.

**LucRo sportic ANKE**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30018	microfibre/ patent nubuck	dark blue	●	36-43	M, W	S, +10, upper



*LucRo sportic*  
**BOOTS**



30062  
Item no.



30068  
Item no.

LucRo sportic **BIANKA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
30062	vegetable nappa/microfibre	grey/dark grey	● ●	36-43	M, W, +10	S, upper
30068	patent nubuck/microfibre	dark blue	●	36-43	M, W	S, +10, upper



**Innovative midsole**

with joint-friendly  
'rebound effect'

**Anti-slip**

TPU outsole

**Four different  
widths**

**Sporty design**

in trendy colours



LucRo  
**sportic**

**FOR  
HIM**



*Peco*

Pages 131–133



*Manuel*

Pages 134–135



*Mario*

Pages 136–137



*Steffen*

Pages 138–139



*Jochen*

Page 141





*LucRo sportic*  
**MEN'S COLLECTION**



**Shoes**  
*with laces/Velcro*

From page 130

**Boots**

From page 140



*LucRo sportic*  
**SHOES**

→ LucRo sportic men's shoes with laces



## LucRo sportic **PECO**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28010	calfskin nappa	black	●	40–50	S, M, W, +10	upper
28011	calfskin nappa/ microfibre 3D mesh	white/grey	○ ●	40–47	M, W	S, +10, upper



**28014**  
Item no.



**Includes spare pair of laces:**  
LucRo flat laces in dark brown,  
140 cm



**28012**  
Item no.

*LucRo sportic* **PECO**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28012	nubuck/carbon look	black	●	40-47	M	S, W, +10, upper
28014	microfibre/ 3D mesh	brown/grey	● ●	40-47	M, W	S, +10, upper



**28019**  
Item no.



✓ **Includes spare pair of laces:**  
LucRo flat laces in grey, 140 cm



**28018**  
Item no.



✓ **Includes spare pair of laces:**  
LucRo flat laces in dark blue,  
140 cm



## LucRo sportic **PECO**

Item no.	Material	Colour	Size	Width ex stock	Delivery time four weeks
28018	microfibre/carbon look	dark blue	● 40-47	M, W	S, +10, upper
28019	microfibre/carbon look	grey	● 40-47	M, W	S, +10, upper

→ LucRo sportic men's shoes with laces



## LucRo sportic **MANUEL**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28020	vegetable nappa/ microfibre	black	●	40-47	M, W	S, +10, upper
28029	calfskin nappa/nubuck	grey/dark blue	● ●	40-47	M, W	S, +10, upper



28028

Item no.




28024

Item no.



KIND VEGETABLE  
LEATHER TANNING



28022

Item no.



## LucRo sportic MANUEL

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28022	calfskin nappa/nubuck	dark blue/grey		40-50	M, W, +10	S, upper
28024	vegetable nappa/ microfibre	brown		40-50	S, M, W, +10	upper
28028	microfibre	dark blue/blue		40-47	M	S, W, +10, upper



## LucRo sportic **MARIO**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28003	calfskin nappa/ 3D mesh	white/blue	<input type="radio"/> <input checked="" type="radio"/>	40-47	M	S, W, +10, upper
28008	microfibre/carbon look	dark blue	<input checked="" type="radio"/>	40-47	M, W	S, +10, upper





28000

Item no.



28001

Item no.



Mario

Model



28009

Item no.

## LucRo sportic MARIO

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28000	calfskin nappa/ 3D mesh	black	●	40-50	S, M, W, +10	upper
28001	calfskin nappa/ 3D mesh	white/grey	○ ●	40-50	M, W, +10	S, upper
28009	calfskin nappa/ honeycomb	grey	●	40-50	M, W, +10	S, upper



## LucRo sportic **STEFFEN**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28030	calfskin nappa/ carbon look	black	●	40-50	S, M, W, +10	upper
28032	calfskin nappa/ patent canvas	dark blue	●	40-50	M, W, +10	S, upper



## LucRo sportic **STEFFEN**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28033	calfskin nappa/ 3D mesh	white/grey	<input type="radio"/> <input checked="" type="radio"/>	40-47	M, W	S, +10, upper
28034	microfibre/ patent canvas	dark brown/brown	<input checked="" type="radio"/> <input type="radio"/>	40-50	M, W, +10	S, upper
28039	microfibre/carbon look	grey	<input checked="" type="radio"/>	40-47	M, W	S, +10, upper



*LucRo sportic*  
**BOOTS**

28042

Item no.

Jochen  
Model



28048

Item no.



## LucRo sportic **JOCHEN**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
28042	microfibre/ canvas patent/nubuck	dark brown/taupe		40–50	M, W, +10	S, upper
28048	microfibre/canvas patent/ carbon look	dark grey/dark blue		40–47	M, W	S, +10, upper



**Anatomically shaped**  
narrow, natural-form last

**Three interchangeable stiffening soles**

**Four different widths**

**Stretch leather**  
(on selected models)

→ LucRo ergonomic women



LucRo  
**ergonomic**

**FOR  
HER**



*Cornelia*

Page 147



*Gisela*

Page 148



*Bärbel*

Page 150



*Helga*

Page 151



*Marlies*

Page 152



*Marianne*

Page 153



*Hannelore*

Page 153



*Carola*

Page 155



*Flora*

Page 157





*LucRo ergonomic*  
**WOMEN'S COLLECTION**



***Sandals***

From page 146

***Shoes***  
*with Velcro/laces*

From page 149

***Boots***

From page 154

***Slippers***

From page 156



*LucRo ergonomic*  
**SANDALS**



LucRo ergonomic **CORNELIA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27420	calfskin nappa/ goatskin nappa metallic	black/anthracite	● ●	36-43	S, M, W	upper
27425	goatskin nappa metallic	sand	●	36-43	S, M, W	upper
27429	goatskin nappa metallic	anthracite/black	● ●	36-43	S, M, W	upper



### LucRo ergonomic **GISELA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27010	calfskin nappa/ embossed patent	black	●	36-43	M, W	S, upper
27011	calfskin nappa	ice	●	36-43	M, W	S, upper
27015	calfskin nappa/nubuck	desert	●	36-43	M	S, W, upper
27018	calfskin nappa/nubuck	dark blue	●	36-43	M	S, W, upper



*LucRo ergonomic*  
**SHOES**

→ LucRo ergonomic women's shoes with Velcro



KIND VEGETABLE  
LEATHER TANNING

## LucRo ergonomic **BÄRBEL**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27042	vegetable nappa	black	●	34-45	S, M, W, +10	upper
27049	antiqued calfskin nappa/ nubuck	grey	●	36-43	M	S, W, +10, upper

27024

Item no.



27028

Item no.



Helga  
Model

27022

Item no.



KIND VEGETABLE  
LEATHER TANNING

## LucRo ergonomic **HELGA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27022	vegetable nappa	black	●	36-43	M, W	S, upper
27024	calfskin nappa/nubuck	cognac	●	36-43	M	S, W, upper
27028	calfskin nappa/nubuck	dark blue	●	36-43	M	S, W, upper

→ LucRo ergonomic women's shoes with Velcro



27061  
Item no.

27060  
Item no.



IDEAL FOR RHEUMATISM



27063  
Item no.

## LucRo ergonomic **MARLIES**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27060	calfskin nappa/ stretch nappa	black	●	36-43	M	S, W
27061	calfskin nappa/ stretch nappa	ice	○	36-43	M	S, W
27063	embossed patent/ stretch nubuck	taupe	●	36-43	M	S, W





IDEAL FOR RHEUMATISM

27070

Item no.

### LucRo ergonomic **MARIANNE**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27070	calfskin nappa/ stretch nubuck	black	●	36-43	M	S, W, +10

→ LucRo ergonomic women's shoes with laces



27031

Item no.

27030

Item no.



### LucRo ergonomic **HANNELORE**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27030	calfskin nappa/nubuck	black	●	36-43	M	S, W, +10, upper
27031	calfskin nappa/ antiqued calfskin nappa	pearl/white	● ○	36-43	M	S, W, +10, upper

*Carola*  
Model

27050  
Item no.



*LucRo ergonomic*  
**BOOTS**

→ LucRo ergonomic women's boots with Velcro



## LucRo ergonomic **CAROLA**

Item no.	Material	Colour	Size	Width ex stock	Delivery time four weeks
27050	calfskin nappa Hydro	black	● 36-43	M, W, +10	S, upper



*LucRo ergonomic*  
**SLIPPERS**



**LucRo Komfort diabetic socks**

for even more cosy comfort at home

**27084**



Item no.



**27089**

Item no.

*LucRo ergonomic* **FLORA**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
27084	felt/ antiqued calfskin nappa	beige		36-43	M	S, W, upper
27089	felt/nubuck	grey/black		36-43	M, W	S, upper



**Anatomically shaped**

narrow, natural-form last

**Three  
interchangeable  
stiffening soles**

**Four different  
widths**

**Stretch leather**

(on selected models)

→ LucRo ergonomic men



LucRo  
**ergonomic**

**FOR  
HIM**





*Lothar*

Page 163



*Jürgen*

Page 165



*Karl*

Page 166



*Reiner*

Page 167



*Heinz*

Page 168



*Manfred*

Page 169



*Detlef*

Page 171



*Ferdinand*

Page 173





*LucRo ergonomic*  
**MEN'S COLLECTION**



***Sandals***

From page 162

***Shoes***  
*with Velcro/laces*

From page 164

***Boots***

From page 170

***Slippers***

From page 172



*LucRo ergonomic*  
**SANDALS**



LucRo ergonomic **LOTHAR**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
26010	calfskin nappa	black	●	40–47	M, W	S, upper
26014	calfskin nappa/nubuck	cognac	●	40–47	M	S, W, upper
26019	antiqued calfskin nappa/ nubuck	grey	●	40–47	M, W	S, upper



*LucRo ergonomic*  
**SHOES**

→ LucRo ergonomic men's shoes with Velcro



26020

Item no.



26024

Item no.



26023

Item no.

Jürgen

Model

## LucRo ergonomic JÜRGEN

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
26020	buckskin	black	●	37-50	M, W, +10	S, upper
26023	buckskin	chestnut	●	37-50	M, W, +10	S, upper
26024	calfskin nappa/nubuck	cognac	●	40-47	M	S, W, +10, upper

→ LucRo ergonomic men's shoes with Velcro



### LucRo ergonomic **KARL**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
26040	calfskin nappa	black	●	37-50	S, M, W, +10	upper
26042	vegetable nappa	black	●	40-47	M	S, W, +10, upper
26043	vegetable nappa	brown	●	40-47	M, W	S, +10, upper



IDEAL FOR RHEUMATISM



LucRo ergonomic **REINER**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
26070	calfskin nappa/ stretch nappa	black	●	37-50	M, W	S, +10
26073	calfskin nappa/ stretch nappa	dark brown	●	40-47	M	S, W, +10

→ LucRo ergonomic men's shoes with laces



## LucRo ergonomic **HEINZ**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
26050	calfskin nappa	black	●	37-50	M	S, W, +10, upper
26053	calfskin nappa	dark brown	●	40-47	M	S, W, +10, upper





26032

Item no.



KIND VEGETABLE  
LEATHER TANNING



26033

Item no.

*Manfred*

Model



KIND VEGETABLE  
LEATHER TANNING

## LucRo ergonomic **MANFRED**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
26032	vegetable nappa	black	●	40-47	M	S, W, +10, upper
26033	vegetable nappa	brown	●	40-47	M	S, W, +10, upper

26060

Item no.

*Detlef*

Model

*LucRo ergonomic*  
**BOOTS**

→ LucRo ergonomic men's boots with Velcro



## LucRo ergonomic **DETLEF**

Item no.	Material	Colour		Size	Width ex stock	Delivery time four weeks
26060	calfskin Hydro/fabric	black	●	40-47	M, W, +10	S, upper



*LucRo ergonomic*  
**SLIPPERS**

→ LucRo ergonomic men's slippers



**LucRo Komfort diabetic socks**

for even more cosy comfort at home

26099  
Item no.

*Ferdinand*  
Model



LucRo ergonomic **FERDINAND**

Item no.	Material	Colour	Size	Width ex stock	Delivery time four weeks
26099	felt/nubuck	grey/black	40-47	M, W	S, upper



*LucRo*  
**ACCESSORIES**



*Shoe care, measuring aid,  
shoehorns, selection guide*

From page 176

*Socks, laces,  
insoles, display unit*

From page 178

# With LucRo READY FOR ANYTHING

## LucRo accessories and collection additions

The comprehensive LucRo shoe range is perfectly complemented with this coordinating range of accessories. Starting with high-quality care products for the ideal protection of leather against moisture and damp, through shoehorns and foot measuring aids, to perfectly fitting specialist socks – LucRo has it all for you.



## LucRo CARE SET

Item no.

230999-000

LucRo care set with contents

### Product features:

- exclusive leather bag with contents, everything also individually available

### Contents:

Two LucRo leather care creams, one polishing cloth, one cream brush, one gloss brush





## LucRo FOOT MEASURING DEVICE

Item no.

210009-001	LucRo foot measuring device
------------	-----------------------------

### Product features:

- up to size 50, for LucRo classic and LucRo kinetic, not suitable for third-party products



## LucRo SHOEHORNS

Item no.	Length	
230575-015	150 mm	no hook
Q86441-016	150 mm	with hook
230575-002	420 mm	with hook
230575-000	500 mm	with hook
230575-010	700 mm	with hook

### Product features:

- metal, chrome-plated



## LucRo ROTATING DISC

Item no.

099952-077	Lucro rotating disc
------------	---------------------

### Product features:

- practical rotating disc for determining correct width



The **Schein Orthopädie Service KG** app can help you find the right LucRo shoe in just a few steps. Download the app at:





## LucRo COMFORT DIABETIC SOCKS

Item no.	Colour		Size	Width ex stock
202100	mottled light grey	●	unisex	35–36 to 47–48
202200	black	●	unisex	35–36 to 47–48
202300	white	○	unisex	35–36 to 47–48
202400	dark brown	●	unisex	35–36 to 47–48

### Product features:

- Comfort socks, semi-plush, hand-knotted
- special padding across the entire sole area
- loose knitted cuff to avoid leg constriction



Item no.	Shoe line	Gender	Starch	Size
229900	LucRo classic	Women	9 mm	34–45
219900	LucRo classic	Men	9 mm	37–50
259900	LucRo kinetic	Women	6 mm	34–45
249900	LucRo kinetic	Men	6 mm	37–50
279900	LucRo ergonomic	Women	6 mm	34–45
269900	LucRo ergonomic	Men	6 mm	37–50
309900	LucRo sportic	Women	6 mm	36–43
289900	LucRo sportic	Men	6 mm	40–50

### Product features:

- Special soft foam with an especially smooth surface to help your foot glide easily into the shoe
- especially suitable for polyneuropathy



Alternatively, Schein also provides NovaPED insert blanks for special indications.  
Find out more at [www.schein.de](http://www.schein.de)



### LucRo LACES

Item no.	Colour	
239477-0	black	●
239477-1	white	○
239477-2	dark brown	●
239477-3	beige	●

**Product features:**

- elastic and available in different lengths (45, 60, 70 and 90 cm)




### LucRo DISPLAY RACK

Item no.	
202020-999	Display rack

**Product features:**

- for LucRo Comfort diabetic socks including 28 pairs of socks (one pair of each colour and size)
- rotates for easy product viewing





'Foot complaints' sources, pages 10–11:

<sup>1</sup>Deutscher Gesundheitsbericht [German health report] 2018, Die Bestandsaufnahme [taking stock], publisher: diabetesDE – Deutsche Diabetes Hilfe [German diabetes aid]

<sup>2</sup>Rheuma in Zahlen – Betroffene Menschen in Deutschland' [rheumatism in numbers – affected people in Germany], publisher: Deutsche Gesellschaft für Rheumatologie e.V., as of November 2015

<sup>3</sup>American Academy of Orthopaedic Surgeons: 'Rheumatoid Arthritis of the Foot and Ankle', at <http://orthoinfo.aaos.org/topic.cfm?topic=a00163> (accessed on 14.2.2017)

<sup>4</sup>ORTHOpress 3/2008, Diagnose & Therapie 'Hallux valgus: die richtige OP für jeden Fuß' [diagnosis and therapy 'hallux valgus: the right op for every foot', pp. 70–71



PEFC-certified

This product is sourced from sustainably managed forests, recycling and verified sources

[www.pefc.de](http://www.pefc.de)



**Schein Orthopädie Service KG**

Hildegardstrasse 5  
42897 Remscheid  
Germany

Tel. +49 (0)2191 910 201  
Fax +49 (0)2191 910 102

lucro@schein.de  
**www.lucro.de**



**Your LucRo dealer:**



099993-007